

## Arthroscopie Systeem Dubbel

REF: 39250

Sterile, single use, non-invasive and non-active, single or double PVC- tubing systems of various length/diameter, with spike(s), tubing clamp(s), connectors, with or without drip chamber, pillow valve or hand pump.

### INTENDED USE:


Intended for washing out a body cavity or surgical sites with a stream of clear sterile fluid.

### INDICATIONS FOR USE

Per the discretion of the qualified medical professional, when washing out a body cavity or surgical site with a stream of clear sterile fluid is warranted.

### INTENDED PATIENT

The intended patient population is restricted to adults, based on the medical condition and other clinical considerations; the qualified medical professional needs to decide if the medical device is suitable for the intended medical procedure.



## PRODUCT

PRODUCT FAMILY:	Irrigation
DURATION OF USE:	Short Term Use: continuous use between 60 minutes and 30 days.
REGULATION:	EU MDD 93/42/EEC
PRODUCT CLASSIFICATION:	Class IIa Rule 2
PRODUCT INFORMATION:	Total Length Tubes: ±280cm Diameter Tubes: Ø4.8x6.8mm Length Silicone: ±8cm Diameter Silicone: Ø5x10mm Open/close

## MATERIALS

MAIN MATERIALS:	LDPE - ABS - Master White 11060 - PVC DEHT 553 - PP - PVC - Silicone
LATEX:	NO
PHthalates, CMR, ENDOCRINE DISRUPTING SUBSTANCES (≥ 0.1% w/w):	NO
PRESENCE OF ANIMAL ORIGIN SUBSTANCES:	YES
SHELF LIFE:	≤ 36 MONTHS

