

Connecting Tube Ch 24, 300 cm

REF: 47750

Single use, sterile, non-active, non-invasive, PVC-tubing of various length/diameter, with or without connectors and/or tubing clamp.

INTENDED USE:

Intended to interface between suction source with collection system and the suction device (probe) during a medical procedure aiming to remove bodily fluids and secretions, surgical tissue debris, or irrigation fluids.

INDICATIONS FOR USE

Per the discretion of the qualified medical professional for any procedure where the surgical site or body cavity has to be cleaned of excess fluid to facilitate observation and/or to clear an obstruction.

INTENDED PATIENT

The intended patient population is unrestricted, based on the medical condition and other clinical considerations; the qualified medical professional needs to decide if the medical device is suitable for the intended medical procedure.



PRODUCT

PRODUCT FAMILY:	Suction Tube
DURATION OF USE:	Short-term: continuous use between 60 minutes and 30 days.
REGULATION:	EU MDD 93/42/EEC
PRODUCT CLASSIFICATION:	Class I Rule 1
PRODUCT INFORMATION:	Length Tube: ±295cm Diameter Tube: Ø5x8mm Length Silicone: ±4cm Diameter Silicone: Ø6x10mm Open/open

MATERIALS

MAIN MATERIALS:	PVC - LN13 - PVC DEHT 903 - Silicone
LATEX:	NO
PHthalates, CMR, ENDOCRINE DISRUPTING SUBSTANCES (≥ 0.1% w/w):	NO
PRESENCE OF ANIMAL ORIGIN SUBSTANCES:	NO
SHELF LIFE:	≤ 36 MONTHS

