

Zuigset Ch 27, 210 cm

REF: 39265

Sterile, single use, surgically invasive, non-active tubing systems. Different types of silicone single lumen tubing of various length/diameters, perforated and partly X-Ray opaque

INTENDED USE:

To be used as a conduit for removal of bodily fluids and secretions, surgical tissue debris, irrigation fluids and gasses, in surgical invasive procedures or through bodily orifices or cavities.

INDICATIONS FOR USE

Per the clinical practitioner's professional discretion for any procedure where the surgical site or body cavity has to be cleaned of excess fluid to facilitate observation and/or to clear an obstruction.

INTENDED PATIENT

The intended patient population is unrestricted, based on the medical condition and other clinical considerations; the qualified medical professional needs to decide if the medical device is suitable for the intended medical procedure.



PRODUCT

PRODUCT FAMILY:	Suction Probe
DURATION OF USE:	Short-term: continuous use between 60 minutes and 30 days.
REGULATION:	EU MDD 93/42/EEC
PRODUCT CLASSIFICATION:	Class IIa Rule 586
PRODUCT INFORMATION:	Length Tube: ±195cm Diameter Tube: Ø6x9mm Open/open

MATERIALS

MAIN MATERIALS:	PVC - LN13 - ABS - PVC DEHT 753
LATEX:	NO
PHthalates, CMR, ENDOCRINE DISRUPTING SUBSTANCES (≥ 0.1% w/w):	NO
PRESENCE OF ANIMAL ORIGIN SUBSTANCES:	NO
SHELF LIFE:	36 MONTHS

