

Zuigslang Ch 30, 200 cm

REF: 48001

Single use, sterile, non-active, non-invasive, PVC-tubing of various length/diameter, with or without connectors and/or tubing clamp.

INTENDED USE:

Intended to interface between suction source with collection system and the suction device (probe) during a medical procedure aiming to remove bodily fluids and secretions, surgical tissue debris, or irrigation fluids.

INDICATIONS FOR USE

Per the discretion of the qualified medical professional for any procedure where the surgical site or body cavity has to be cleaned of excess fluid to facilitate observation and/or to clear an obstruction.

INTENDED PATIENT

The intended patient population is unrestricted, based on the medical condition and other clinical considerations; the qualified medical professional needs to decide if the medical device is suitable for the intended medical procedure.



PRODUCT

PRODUCT FAMILY:	Suction Tube
DURATION OF USE:	Short-term: continuous use between 60 minutes and 30 days.
REGULATION:	EU MDD 93/42/EEC
PRODUCT CLASSIFICATION:	Class I Rule 1
PRODUCT INFORMATION:	Length Tube: ±200cm Diameter Tube: Ø7x10mm Open/open

MATERIALS

MAIN MATERIALS:	PVC - LN13 - PVC DEHT 603
LATEX:	NO
PHthalates, CMR, ENDOCRINE DISRUPTING SUBSTANCES (≥ 0.1% w/w):	NO
PRESENCE OF ANIMAL ORIGIN SUBSTANCES:	NO
SHELF LIFE:	36 MONTHS

