

## Suction Set Ch 24, 200 cm

REF: 59007

Sterile, single use, surgically invasive, non-active tubing systems. Different types of silicone single lumen tubing of various length/diameters, perforated and partly X-Ray opaque

### INTENDED USE:


To be used as a conduit for removal of bodily fluids and secretions, surgical tissue debris, irrigation fluids and gasses, in surgical invasive procedures or through bodily orifices or cavities.

### INDICATIONS FOR USE

Per the clinical practitioner's professional discretion for any procedure where the surgical site or body cavity has to be cleaned of excess fluid to facilitate observation and/or to clear an obstruction.

### INTENDED PATIENT

The intended patient population is unrestricted, based on the medical condition and other clinical considerations; the qualified medical professional needs to decide if the medical device is suitable for the intended medical procedure.



## PRODUCT

PRODUCT FAMILY:	Suction Probe
DURATION OF USE:	Short-term: continuous use between 60 minutes and 30 days.
REGULATION:	EU MDD 93/42/EEC
PRODUCT CLASSIFICATION:	Class IIa Rule 586
PRODUCT INFORMATION:	Length Tube: ±195cm Diameter Tube: Ø5x8mm Open/open

## MATERIALS

MAIN MATERIALS:	PETG SKYGREEN S2008 - PVC - LN13 - PVC DEHT 903
LATEX:	NO
PHthalates, CMR, ENDOCRINE DISRUPTING SUBSTANCES (≥ 0.1% w/w):	NO
PRESENCE OF ANIMAL ORIGIN SUBSTANCES:	NO
SHELF LIFE:	36 MONTHS

